




Active Magnetic Field Canceller

AMC-331

Protects Electron Beam application equipment from environmental magnetic interferences

High performance Active Magnetic Field Canceller **AMC-331**



- Scanning Electron Microscope (SEM)
- Transmission Electron Microscope (TEM)
- Focused Ion Beam (FIB)
- Wafer Inspection System (CD-SEM)
- Electron Beam Lithography (EBL)
- Mass Spectrometer (MS)
- Others , EB application equipment



Amplifier & Controller
(DM-70)

The high performance active magnetic field canceller **AMC-331**

Magnetic field fluctuations due to indoor wiring, cars, trains and elevators adversely affect the electron beam as a disturbance of magnetic field, which is an obstacle to using these devices. The tool that solves this problem is the active magnetic field canceller system. The active magnetic field canceller system "AMC-331" is composed of active magnetic field canceller "DM-70", with the magnetic field sensor and the control coil. It helps to suppress the adverse effect of the magnetic field disturbance and improve the image of the electron microscope and the processing accuracy of the electron beam machining apparatus.


Basic Theory(Helmholtz coil)

■ Causes of Magnetic Field Fluctuations

- Movement of magnetic bodies (automobiles, elevators, trains, etc.)
- Electric power supply noise (50 or 60Hz)
- Peripheral equipment

■ Helmholtz Coil

Helmholtz coils are a pair of parallel coils of the same number of turns connected in series. They will generate the magnetic field along the axis through the center of the two coils when the electric current flows in the coils.




To generate compensation magnetic field to cancel the magnetic disturbance.

System configuration

■ Configuration


- ① Amplifier & Controller
- ② Magnetometer
- ③ Control Software
- ④ Magnetic field generating coil

■ System block diagram




■ Magnetic field measurement data

Contour map of magnetic flux density



Parallel magnetic field from Helmholtz coil



※PC for monitor is not included in the product.
 ※The cage (coil cover) varies according to the object and the demand of customers. It is optional, and should be ordered separately corresponding to the application.

Features

- High Performance
 - High damping performance of maximum 40~60dB
- Low Cost
 - Cost of 1/10 against passive shield solution
- Easy Installation
 - Easy installation just for arrangement of coils around the equipment to be protected
- Optimum Design
 - Coil design in accordance with the target tool.
- Digital Control
 - Pin-point control of several magnetic field disturbances based upon the frequency analysis.

Control Procedure

Quasi DC + Typical AC component control + Independent feedback control for each axis

Control Range

$\pm 15\mu\text{T}$

Control Frequency Range

DC~1kHz

Maximum Reduction Ratio


40~60dB

Application


- Scanning Electron Microscope (SEM)
- Transmission Electron Microscope (TEM)
- Focused Ion Beam(FIB)
- Wafer Inspection System (CD-SEM)
- Electron Beam Lithography (EBL)
- Mass Spectrometer(MS)

Experimental results


1) Experimental results at 200Hz disturbance




2) Experimental results with environmental disturbance



SEM image improvement by AMC-331 OFF / ON



Coil Cable Installation Option



In addition, we can provide installation proposals and dual sensor solution according to specifications and requirements. In summary, AMC-331 protects the equipment from electro magnetic interferences.

AMC-331 Specifications

Description of Whole System	
Model	AMC-331
Main Configuration	Controller
	Magnetometer
	Control Software
	Magnetic field generating coil
Control Method	Three axes (X,Y,Z) Simultaneous Control
	Quasi DC+Typical AC component control
	Independent Feedback Control for each axis
Control Range	$\pm 15\mu\text{T}$ (at 2 x 2 x 2m)
Control Frequency Range	Quasi DC 1kHz
Attenuation Factor	Maximum 40~60dB

Controller	
Model	Digital Controller DM-70
Main Configuration	Control / Processor(CPU-DSP)Board, Sampling(AD-DA)Board
Sampling Frequency	2kHz
A/D and D/A resolution	A/D 18bit, D/A 16bit
Input Channel	3ch / 3 axis of Fluxgate Sensor (Dual Sensor for Option)
Output Channel	3ch / 3 axis of output coils, BNC for monitoring (3ch / 3 axis, Dual Sensor for Option)
Display	Digital display for X/Y/Z EMI value (At sensor point, AC/DC switchable, μT_{p-p})
AC Power	AC 100-240V 50/60Hz 300W
Dimension	430(W) x 460(D) x 44.5(H)mm (Except protruding parts)
Weight	4.5kg

Magnetometer	
Sensor Type	Fluxgate
Frequency Range	DC ~ 1kHz
Measurement range	$\pm 100\mu\text{T}$ (1G)

Control Software	
Model	AMC Mon 1.2
Main Function	Real-time Wave Monitor(Time Series of FFT)
	Control Parameters Adjusting
	Wave Data Recording to a Log File

Coils for Generation Magnetic Field	
Coil Structure	3 Axes Helmholtz-Coil Cage
Output Coil Current	$\pm 2\text{A}$ /ch Max

	The caution for safety	Please read the Instruction Manuals carefully before use.
--	-------------------------------	---

This catalogue is for the product as of Oct.2020. The specification may be modified without notice.



Vibration Control Technology
特許機器株式会社

TOKKYOKIKI

Head Office 10-133, Minami-Hatsushima-Cho, Amagasaki City, Hyogo, 660-0833, Japan
(Factory) Tel: (+81) 6 - 6487 - 3939 Fax: (+81) 6 - 6487 - 3947

Tokyo Office 2-5-15, Higashi-Kanda, Chiyoda-ku, Tokyo, 101-0031, Japan
Tel: (+81) 3 - 6831 - 0011 Fax: (+81) 3 - 6831 - 0008

精密振動系商品のご相談は、豊富な経験・実績で信頼できる特許機器へ [https:// www.tokkyokiki.co.jp](https://www.tokkyokiki.co.jp)

Cat No.20.10.05EN